



**Statistical Analysis Plan for Masimo Study titled
“Calibration and Validation of Masimo’s O3
Regional Oximetry Device in Neonates, Infants
and Children Undergoing Cardiac
Catherization”**

Protocol Date: 9 MAY 2018

NCT Number: NCT03123354

Statistical Analysis Plan for Masimo Study titled “Calibration and Validation of Masimo’s O3 Regional Oximetry Device in Neonates, Infants and Children Undergoing Cardiac Catherization”

Data Analysis

We compute error statistics, including mean bias, precision, standard deviation of bias, and absolute root-mean-squared-error accuracy A_{rms} , accounting for within-subject correlations.

Applicable Analysis Methods

Measurement error or bias for the k th sample $e(k)$ is defined as

$$e(k) = \text{rSO}_2(k) - \text{SavO}_2(k).$$

Mean bias over N data points is defined as

$$\mu = \frac{1}{N} \sum_{k=1}^N e(k)$$

The A_{rms} can be calculated as

$$A_{\text{rms}} \approx \sqrt{\mu^2 + s^2}$$

where s is the estimated standard deviation of bias, accounting for within-subject correlation, which is calculated using the Bland and Altman method for repeated measures. The 95% confidence interval of mean bias was calculated as $(\mu - 2s/\sqrt{N}, \mu + 2s/\sqrt{N})$.