

TITLE: Bone healing at non-submerged implants installed with different insertion torques. A split mouth histomorphometric randomized controlled trial.

AUTHORS: Yoshiyuki Amari, DDS*, Adriano Piattelli, MD, DDS, DrHC, DrHC**, Karol Ali Apaza Alccayhuaman DDS***, Natalia Fortich Mesa, DDS, MSc§, Mauro Ferri, DDS#, Giovanna Iezzi, DDS, PhD##, Daniele Botticelli, BMBD, PhD°

* Research affiliate at ARDEC Academy, Ariminum Odontologica, Rimini, Italy;

** Full Professor, Department of Medical, Oral and Biotechnological Sciences, University of Chieti-Pescara, Chieti, Italy; Chair of Biomaterials Engineering, Catholic University of San Antonio of Murcia (UCAM), Murcia, Spain;

*** Research affiliate at ARDEC Academy, Ariminum Odontologica, Rimini, Italy;

§ Docente con Función de Coordinación en Investigación, Programa de Odontología, Corporación Universitaria Rafael Nuñez, Cartagena de Indias, Colombia;

Docente Auxiliar II, Facultad ciencias de la salud, Programa de Odontología, Universidad Rafael Núñez, Cartagena de Indias, Colombia;

Associate Professor, Department of Medical, Oral and Biotechnological Sciences, University of Chieti-Pescara, Chieti, Italy;

° Private Practice, Principal Researcher at ARDEC Academy, Ariminum Odontologica, Rimini, Italy

Corresponding author:

Karol Ali Apaza Alccayhuaman

Viale Giovanni Pascoli 67, 47923 Rimini (RN), Italy

Tel.: +39-0541-393444, Fax: +39-0541-397044, E-mail: caroline7_k@hotmail.com

Cartagena 20th April 2017

ABSTRACT

Objectives: to evaluate histomorphometrically the healing at implants installed with standard or very low insertion torque values.

Material and methods: Twelve volunteer patients will be recruited, and two screw-shaped titanium devices will be installed in the distal segments of the mandible using insertion torque values of either <10 Ncm or ~30 Ncm. The implants will be left to heal in a non-submerged fashion. After 8 weeks, biopsies will be retrieved and ground sections will be prepared for histological evaluation.

INTRODUCTION

Very high insertion torque values have been recommended when immediate loading is applied to implants.¹ However, it has been shown that similar clinical outcomes may be achieved even with insertion torque values ≤ 15 Ncm when implants are splinted together.² Nevertheless, in a consensus conference, insertion torque values comprised between 20 and 45 Ncm were recommended for implants immediately loaded with single crowns.³ Moreover, in animal experiments, the highest rate of osseointegration was observed when insertion torque of 30-35 Ncm were applied.⁴⁻⁶

Different insertion torque values have been tested both in human⁷⁻⁹ and in animals studies.^{1,4-6,10-12} A systematic review with meta-analyses summarized data from both animal and human studies and two groups were identified based on the insertion torque values: >25 Ncm or <30 Ncm.¹³ No differences were found between the two groups in terms of implant survival rate or marginal bone loss. Nevertheless, in a randomized clinical study,⁷ after 12 months of healing, implants installed with torque ≥ 50 Ncm presented higher marginal bone loss compared to implants with torque included between 20-50 Ncm. Very low insertion torques have been also tested both in animals and humans. After 4 months of healing, in an experimental study in dogs,⁵ an osseointegration >55% was

obtained at implant installed with torque close to 0 Ncm. In a clinical study,⁸ eleven implants were installed with a torque <10 Ncm. After 4 to 6 months, a reversal torque of 35 Ncm was applied and all implants showed a good stability.

However, it has to be considered that implants presenting a rotational instability during implant installation resulted in a lower survival rate compared to implants installed with higher insertion torque.^{14,15}

Due to the contradictory outcomes on the influence of the torque on osseointegration and a lack of histological data in humans, there is a need of more evidences that may support the clinicians in the decision making when an unintentional low insertion torque occur at implants during the daily practice.

Hence, the aim of the present study was to evaluate histomorphometrically the healing at implants installed with standard or very low insertion torque values.

MATERIALS & METHODS

Patient selection

The Declaration of Helsinki on medical protocols and ethics will be followed. The protocol has the approval of the Ethical Committee of the Corporación Universitaria Rafael Núñez, Cartagena de Indias, Colombia with protocol #04-2014 on October 8, 2014. All treatments will be performed in that institution starting from November 2016. All surgical procedures and possible complication will be clearly explained to each participant and a written informed consent was signed by all patients. The study will be reported according to the CONSORT guidelines.

For sample calculation, data from a dog experiment were used⁵ and a sample of 12 subjects was calculated for matched pairs considering a difference in bone-to-implant contact of 10% being clinically relevant. A power of 0.8 and $\alpha=0.05$ were used.

The following inclusion criteria will be required:

- (i) presence of at least two edentulous zone in the posterior segment of the mandible
- (ii) ≥ 25 years of age
- (iii) smoking ≤ 10 cigarettes per day
- (iv) Good general health
- (v) No contraindication for oral surgical procedures.
- (vi) Not being pregnant.

The following exclusion criteria will be adopted:

- (i) Presence of systemic disorders
- (ii) Chemotherapy or radiotherapy;
- (iii) Smokers >10 cigarettes per day
- (iv) Previous bone augmentation procedures in the region.

Device

Customized solid titanium screw-shaped devices will be used (Sweden & Martina, Due Carrare, Padua, Italy). The devices have an intraosseous portion with a moderately rough surface¹⁶ (ZirTi® surface, Sweden & Martina, Due Carrare, Padua, Italy). The intraosseous portions will be 4 mm long, with a diameter of 2.65 mm at the apical aspect and 3 mm at the coronal margin. A polished neck 2.4 mm long junction represents the transmucosal portion of the implant.

Randomization

Each patient will receive two mini-implants, installed in recipient sites prepared either to reach torque of <10 Ncm or ~ 30 Ncm. The two recipient sites will be selected prior the surgery, while the type of site preparation will be randomly decided. The randomization will

be performed electronically (randomization.com) by a researcher neither involved in the selection of the patients nor in installation of the devices (DB). Sealed opaque envelopes will be prepared and opened at the time of surgery by an author not involved in the surgery (DB). Surgeon and patients will be blinded.

The histological slides will be coded so that the assessor will be blinded about the type of site preparation when the histological slides will be evaluated.

Clinical procedures

The surgical procedures will be performed by an expert surgeon. After the injection of local anesthesia, small crestal and releasing incisions will be performed in the distal segments of the mandible, and small full-thickness muco-periosteal flaps will be elevated. The test sites (<10 Ncm) will be over-prepared with drills of larger diameter compared to those used at the standard sites (~30 Ncm). All site preparations will be performed deeper compared to the length of the implant so that the implant apex cannot reach the bottom of the osteotomy. This, in turn, means that the final torque will be produced by the lateral pressure against the bone walls, as described in a previous animal experiment.⁴ The mini-implants will be subsequently installed and the final insertion torque will be measured with wrench calibrated on a newton-meter (Newton; Leader Italia, Cinisello Balsamo, MI, Italy). A cover screw will be placed on the top of the mini-implants, and the flaps will be sutured allowing a non-submerged healing.

Maintenance

Antibiotics (amoxicillin 875 mg/clavulanic acid 125 mg twice a day for 6 days) and non-steroidal anti-inflammatory drugs as needed (Ibuprofen 400 mg) will be prescribed. Mouth rinses with 0.12% chlorhexidine three times a day for 10 days will be also recommended.

After 7 days, the sutures will be removed, and the patients will be recalled every 2 weeks after surgery.

Biopsies

After 8 weeks of healing the patients will be recalled to the clinic for biopsies retrieval. Full-thickness flaps will be elevated and biopsies including the mini-implants will be retrieved using the trephine in an eccentric position to reduce the size of the donor site and obtain sufficient hard tissue at least at one side of the biopsies.¹⁷

Histological preparation of the biopsies

The biopsies will be washed in saline solution and immediately stored in 10% buffered formalin. The histological process will be performed in the Histology Laboratory for Hard Tissues at the University of Chieti-Pescara, Italy. All biopsies will be first dehydrated in an ascending series of alcohol and subsequently embedded in resin (Technovit® 7200 VLC; Kulzer, Wehrheim, Germany). After polymerization, the biopsies will be sectioned following the longitudinal axis using a precision diamond disk. Specimens of ~150 microns of width will be obtained that will be afterwards ground to ~30 µm of width. A staining with acid fuchsine and toluidine blue will be applied.

Histomorphometric evaluation

The histomorphometric evaluations will be performed twice by an expert examiner, and mean values will be used. The examiner will be blinded, and no indications will be reported on the histological slides that may allow the identification of groups. All histological analyses will be performed in ARDEC facilities (Ariminum Odontologica, Rimini, Italy) using an Eclipse Ci microscope (Nikon Corporation, Tokyo, Japan) that will be coupled with a digital video camera (Digital Sight DS-2Mv, Nikon Corporation, Tokyo, Japan). The

measurements will be taken using the software NIS-Elements D 4.10 (Laboratory Imaging, Nikon Corporation, Tokyo, Japan). The percentages of new bone, pre-existing (old) bone, soft tissue (marrow spaces, Haversian canals, BMUs canals), bone debris/ clot remnants, and marrow spaces (soft tissue) will be evaluated both in contact with the implant surface and surrounding the implant to a distance of 0.4 mm from the surface. All measurements will be performed at x200 magnification from the most coronal contact of bone to the implant surface (B) to the apical extension of the tissues (A). The total mineralized bone will be assessed as sum of new and old bone. For the morphometrical measurements of the tissues surrounding the implant surface (tissue density), a point counting procedure will be adopted, using a lattice with squares of 50 microns superimposed over the histological image, at a magnification of x200.

Outcomes measures

Primary outcome: New bone in contact with the implant surface

Description: The percentages of new bone will be evaluated both in contact with the implant surface

Timeframe: After 8 weeks of healing

Secondary outcome: Bone density around the implant surface

Description: The total mineralized bone will be assessed as sum of new and old bone

Timeframe: After 8 weeks of healing

Other pre-specified outcomes: pre-existing (old) bone, soft tissue (marrow spaces, Haversian canals, BMUs canals), bone debris/ clot remnants, and marrow spaces (soft tissue)

Description: The percentages of pre-existing (old) bone, soft tissue (marrow spaces, Haversian canals, BMUs canals), bone debris/ clot remnants, and marrow spaces (soft tissue) will be evaluated all in contact with the implant surface

Timeframe: After 8 weeks of healing.

Data analysis

The primary variable will be new bone in contact with the implant surface. New bone (bone density) around the implant surface will be the secondary variable.

Mean values and standard deviations as well as 25th, 50th (median), and 75th percentiles will be calculated for each outcome variable. Means, standard deviations and 95% confidence intervals of the differences between test and control sites will be calculated for each variable analyzed.

A Wilcoxon test will be used to analyze differences between sonic and drill groups. The level of significance will be set at α 0.05.

ACKNOWLEDGEMENTS

This study will be economically supported by Sweden & Martina, Due Carrare, Padua, Italy, by ARDEC Academy, Ariminum Odontologica Srl, Rimini, Italia.

DISCLOSURE: The Authors claim to have no financial interest, either directly or indirectly, in the products or information listed in the article.

REFERENCES

1. Trisi P, Todisco M, Consolo U, Travaglini D. High versus low implant insertion torque: a histologic, histomorphometric, and biomechanical study in the sheep mandible. *Int J Oral Maxillofac Implants*. 2011;26:837-849.
2. Cesaretti G, Botticelli D, Renzi A, et al. Radiographic evaluation of immediately loaded implants supporting 2-3 units fixed bridges in the posterior maxilla: a 3-year follow-up prospective randomized controlled multicenter clinical study. *Clin Oral Implants Res*. 2016;27:399-405.
3. Gallucci GO, Benic GI, Eckert SE, et al. Consensus statements and clinical recommendations for implant loading protocols. *Int J Oral Maxillofac Implants*. 2014;29 Suppl:287-290.
4. Pantani F, Botticelli D, Garcia IR Jr, et al. Influence of lateral pressure to the implant bed on osseointegration: an experimental study in dogs. *Clin Oral Implants Res*. 2010;21:1264-1270.
5. Rea M, Lang NP, Ricci S, et al. Healing of implants installed in over- or under-prepared sites--an experimental study in dogs. *Clin Oral Implants Res*. 2015;26:442-446.
6. Rea M, Botticelli D, Ricci S, et al. Influence of immediate loading on healing of implants installed with different insertion torques--an experimental study in dogs. *Clin Oral Implants Res*. 2015;26:90-95.
7. Barone A, Alfonsi F, Derchi G, et al. The effect of insertion torque on the clinical outcome of single implants: a randomized clinical trial. *Clin Implant Dent Relat Res*. 2016;18:588-600.
8. Verardi S, Swoboda J, Rebaudi F, Rebaudi A. Osteointegration of tissue-level implants with very low insertion torque in soft bone: a clinical study on SLA surface treatment. *Implant Dent*. 2018;27:5-9.
9. Verrastro Neto A, Andrade R, Corrêa MG, ET al. The impact of different torques for the insertion of immediately loaded implants on the peri-implant levels of angiogenesis- and bone-related markers. *Int J Oral Maxillofac Surg*. 2018;47:651-657.
10. Duyck J, Corpas L, Vermeiren S, et al. Histological, histomorphometrical, and radiological evaluation of an experimental implant design with a high insertion torque. *Clin Oral Implants Res*. 2010;21:877-884.
11. Duyck J, Roesems R, Cardoso MV, et al. Effect of insertion torque on titanium implant osseointegration: an animal experimental study. *Clin Oral Implants Res*. 2015;26:191-196.
12. Mukhtadar AK, Gangaiah M, Chrcanovic BR, Chowdhary R. Evaluation of the effect of self-cutting and nonself-cutting thread designed implant with different thread depth on variable insertion torques: An histomorphometric analysis in rabbits. *Clin Implant Dent Relat Res*. 2018;20:507-514.
13. Berardini M, Trisi P, Sinjari B, et al. The effects of high insertion torque versus low insertion torque on marginal bone resorption and implant failure rates: a systematic review with meta-analyses. *Implant Dent*. 2016;25:532-540.

14. Balshi SF, Wolfinger GJ, Balshi TJ. A retrospective analysis of 44 implants with no rotational primary stability used for fixed prosthesis anchorage. *Int J Oral Maxillofac Implants*. 2007;22:467-471.
15. Rodrigo D, Aracil L, Martin C, Sanz M. Diagnosis of implant stability and its impact on implant survival: a prospective case series study. *Clin Oral Implants Res*. 2010;21:255-261.
16. Caneva M, Lang NP, Calvo Guirado JL, et al. Bone healing at bicortically installed implants with different surface configurations. An experimental study in rabbits. *Clin Oral Implants Res*. 2015;26:293-299.
17. Ferri M, Lang NP, Angarita Alfonso EE, et al. Use of sonic instruments for implant biopsy retrieval. *Clin Oral Implants Res*. 2015;26:1237-1243.

